



wave

DESIGNED TO REDUCE THE
APPEARANCE OF CELLULITE &
ENHANCE BODY CONTOURING

BODY CONTOURING
CELLULITE REDUCTION
SKIN TIGHTENING
SKIN ELASTICITY
COLLAGEN REGENERATION



THE PERFECT COMBINATION FOR SKIN IMPROVEMENT DIRECTLY FOLLOWING FAT REDUCTION TREATMENTS.

The Z Wave is a clinically proven treatment for cellulite and body contouring. The treatment targets and rejuvenates the skin and underlying soft tissues to smooth out cellulite, tighten loose skin, and reduce the appearance of stretch marks

The treatment works by breaking down fat under the skin and stimulating lymphatic drainage to remove excess fat materials and toxins from the body. It also aids in collagen production for skin that appears smoother, healthier, and more youthful.

INDICATIONS

SKIN LIFTING & TIGHTENING

- Skin Laxity
- Skin Tightening
- Face & Body Contouring

SKIN CONTOURING

- Uneven Skin Tone & Texture
- Sagging Skin & Laxity

THE Z-WAVE PROCEDURE

CONSULTATION

- Your practitioner will advise on suitability and your treatment outcome objectives.
- Treatment zone markings will be made
- For body treatments, measurements and treatment area pictures will be taken.

DURING THE TREATMENT

- Ultrasound gel is applied to the area.
- The handpiece will then be lightly applied and the Z Wave treatment will begin.

AFTER THE TREATMENT

- The skin may appear red in the treated area but this usually goes within 24 hours.

RESULTS

- A course of 6-8 body treatments is recommended 1-2 times per week
- Improvement in the skin's texture is seen after just four treatments. However, the full results will not be seen until the course has been completed. Further improvements continue 3 months post treatment as collagen production and circulation improves.

The Z-Wave treatment areas:

ABDOMEN

LEGS

UPPER ARMS

BUTTOCKS



QUICK TREATMENT



NON-INVASIVE TREATMENT



NO DOWNTIME



GENTLE & PAINLESS



FREQUENTLY ASKED QUESTIONS?

Who is the Z Wave for?

The Z Wave is recommended for patients who are experiencing resistant cellulite, especially on common problem areas such as the thighs, arms, and buttocks. It is also ideal for patients who have loose, saggy, or crepey skin, and wish to regain smoother, healthier skin.

Is the treatment painful?

The treatment is extremely comfortable.

How long does each session last?

Each treatment area would take approximately 8 minutes.

Are the treatments long-lasting?


Regular refresher sessions, exercise and a balanced nutrition, will ensure long-lasting results.

What are the Contraindications?


Implantable metal electrical device such as a pacemaker, cardioverter or defibrillator, Pregnancy or breastfeeding, Keloids and history of skin cancer in the areas to be treated.

HOW DOES IT WORK?


The Z Wave unit transmits powerful sound waves onto and through the skin. The energy then spreads out radially. The body responds to the Z Wave sound waves with increased metabolic activity in the form of lymphatic drainage and increased blood circulation.



The high energy radial shockwaves consist of two different parts – positive pressure pulse and a comparatively small tensile wave component. The shockwave squeezes the surrounding tissue and breaks the grid structure of the molecules.



The tensile wave leads to a collapse of the gas bubbles within the fat structures (cavitation), which re-expand after the high energy radial shockwaves and thus ensure a destabilization of the fat structures.



POWERED BY



NOVUSMEDICAL



0333 939 7933

INFO@NOVUSMEDICALUK.COM

WWW.NOVUSMEDICALUK.COM

